

# *insta***BAG**<sup>®</sup>

Filter blender bag with dehydrated media

## Your diluent ready in minutes!



# interscience

# Ready-to-use media for easy preparation

*instaBAG*<sup>®</sup> is a filter bag which includes a dehydrated pre-dosed media, allowing a 2-in-1 process. Skip the time-consuming media preparation.



1 Place the sample in an *instaBAG*<sup>®</sup>



2 Add sterile and deionized water



3 Blend the *instaBAG*<sup>®</sup> with the BagMixer<sup>®</sup>

➔ The sample is ready to be pipetted

## All the advantages of a filter bag

Debris are blocked by the filter during blending

## Long shelf life

Easy pouch opening

Zip pouch with desiccant: optimal conservation after opening

Light and humidity proof aluminum pouch



## Just add water

Ready-to-use bag with pre-dosed dehydrated media

Complete dissolution during blending

Non-woven filter identical to the bag filter

Available in BPW 90 mL and BPW 225 mL



PATENTS PENDING

SHELF LIFE 12 MONTHS

STORAGE 6°C to 30°C

## Space saving




100  
*instaBAG*



100  
glass bottles

## No waste

100 %  
used  
3 *instaBAG*  
for 3 analyses.



13 %  
used  
1 (5L) pouch  
for 3 analyses.  
The leftover liquid  
is thrown away.

## Less garbage



*instaBAG*  
+ diluent + sample



Media bottle  
+ bag + diluent + sample

## Easy preparation



with *instaBAG*.



with media preparation

## Green

*instaBAG* : no water transported



Diluent pouches




## Long shelf life

*instaBAG*



Media bottle (BPW)



# Technical specifications

	<i>instaBAG</i> ® BPW 90	<i>instaBAG</i> ® BPW 225 mL
Reference	114 090	114 225
Type of bag	Bag with lateral filter	
Composition of the bag	MultiLayer®: reinforced multicoated complex	
Type of filter	Non-woven	
Filter porosity	< 250 microns	
Type of dehydrated media	Buffered peptone water	
Weight of dehydrated media	1.8 g	4.5 g
Sample weight	10 g (± 5 %)	25 g (± 5 %)
Volume of water to be added (sterile and deionized)	90 mL (± 2 %)	225 mL (± 2 %)
Blending time	1 mn	
Bag dimensions	190 x 300 mm	
Aluminum pouch of	10 bags	
Shelf life	12 months	
Storage conditions	Temperature from + 6°C to + 30°C	
Treatment	Gamma ray (10 to 25 KGy)	
Standards	ISO 11133 - ISO 11290-2 - ISO 215281 - ISO 7218 - ISO 6579 - ISO 6887- FDA-BAM	
Box of	100 bags	

MSDS available on request

Ask for free samples

**Certified production**

Products made for INTERSCIENCE by Interlab, an ISO 9001 certified company.



## The winning combination



**DiluFlow**® Gravimetric dilutors



**instaBAG**®



**BagMixer**® Lab blenders

Your local distributor

**interscience**

30, chemin du Bois des Arpents - 78860 St Nom - FRANCE  
T: +33 (0)1 34 62 62 61 - Email: info@interscience.com

**interscience** USA & CANADA

32 Cummings Park - Woburn, MA 01801 - USA  
P: +1 781 937 0007 - F: +1 781 937 0017 - Email: sales.usa@interscience.com

**interscience** CHINA

上海市徐汇区肇嘉浜路798号坤阳大厦1903室 - 200030  
电话: +86 (0)21-64739390 - +86 189 3097 0733 - 邮址: sales.china@interscience.com

**interscience** SOUTH-EAST ASIA

The Centropod, 80 Changi Rd - 05-07 Singapore 419715  
T: +65 6909 0825 - M: +65 8118 5178 - E-mail: sales.asia@interscience.com

[www.interscience.com](http://www.interscience.com)