

# Ferrous Sulphate Fertilizer

## DIGESTION OF A FERROUS SULPHATE FERTILIZER USING COLDBLOCK™ DIGESTION TECHNOLOGY

*Author: Lorie-Anne Doig*

### Introduction

This application note will focus on the digestion of Ferrous Sulphate Fertilizer pellets using ColdBlock™ Digestion CB12L Technology.

### Method

A sample of Ferrous Sulphate Fertilizer pellets was digested using the following method:

- 0.5g of sample was weighed and placed into a ColdBlock™ Digestion vessel.
- 20ml of 1:1 HNO<sub>3</sub>:HCl was added.
- Sample was digested at 100% power for 10min.
- The chiller was set to -5°C.
- Sample was allowed to cool.
- Sample was normalized to 50mL using 2% HNO<sub>3</sub>.

### Instrument

ColdBlock™ Digestion CB12L Technology.

### General

This procedure is specific for the sample digested and may need modification for different samples to achieve the desired result.